

INSTRUCTIONS FOR USE TUBES FOR OXYGEN THERAPY

Tubes for oxygen therapy MACLYSA.

Practical guidance for placement and use.

Caution:

Use only by qualified personnel. Read these instructions carefully before using the product.

Secure system connections.

Purpose:

The tubes for oxygen therapy MACLYSA have been designed to allow the flow of oxygen from an oxygen therapy equipment up to a certain distance (the length of the tubes) to treat or relieve the patient's breathing problems.

They are designed to ensure that oxygen continues to flow even when the tube is kinked or compressed.

Description:

The tubes for oxygen therapy MACLYSA is made up of three pieces: a PVC tube and two universal connectors attached to each end of the tube de PVC. One of the connectors is connected to the nipple of the oxygen therapy machine and the other to the mask or nasal goggle that will be placed in the patient's nostrils.

Use of tubes for oxygen therapy:

Exclusive use to conduct the flow of oxygen from an oxygen therapy device to the patient.

The tubes for oxygen therapy MACLYSA will be manufactured, packaged and consumed in a NON-STERILE manner.

Contraindications:

Hypersensitivity to any of the constituent elements of the tube.

Any alteration of the properties of the product or of its performance must make it desist from its use, and the anomaly must be immediately communicated to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Warnings:

Single patient use

Store under normal conditions for medical supplies.

Store without exceeding 40°C.

Use a cloth and water for cleaning

Do not reprocess

Dispose of the product as plastic material in accordance with current legislation

Manufactured by:

MACLYSA, S.L.

Polígono Industrial Las Viñas C/ Murcia, 69

13420 Malagón (Ciudad Real) Spain

Phone number: (+34) 926 80 21 78

e-mail: info@maclysa.com



**SYMBOLS LABELS
TUBES FOR OXYGEN THERAPY MACLYSA**

	Batch code
	Use-by date
	Manufacturer
	Catalogue number
	CE-marking in accordance with Regulation (EU) 2017/745
	Non-sterile
	Consult instruction for use
	PVC product
	Upper limit of temperature
	Single patient multiple use
	Medical device

Manufactured by:

MACLYSA, S.L.

Polígono Industrial Las Viñas C/ Murcia, 69

13420 Malagón (Ciudad Real) Spain

Phone number: (+34) 926 80 21 78

e-mail: info@maclysa-sl.com

